

Glyoxylic acid monohydrate test method

色谱柱型号:NG-SA-4A,250*4.6mm。淋洗液:4.8mol/L

Na₂CO₃. 3.6mol/L NaHCO₃

色谱柱柱温:36摄氏度

使用电导检测器检测

Thanslate into English:

Chromatographic column type: NG-SA-4A,250*4.6mm。

Leacheate:4.8mol/L Na₂CO₃. 3.6mol/L NaHCO₃

Chromatographic column temperature: 36°C

Using a conductivity detector to test.